

HemosIL®

Providing the Right Solutions for Your APTT Testing

Activated Partial Thromboplastin Time

*Sensitivities for*  
**Heparin, Intrinsic Factors**  
**and Lupus Anticoagulants**  
*To Meet Your Lab's Requirements*

# Providing the Right Solutions for Your APTT Testing

## Introduction

The Activated Partial Thromboplastin Time (APTT) is an important assay for screening deficiencies and inhibitors of the intrinsic pathway of the coagulation cascade. In an APTT reagent, negatively charged particles (e.g., silica or Ellagic acid) are mixed with specific phospholipids and buffers to make an environment within which intrinsic plasma proteins are easily activated.

After plasma and the APTT reagent incubate together, Calcium is added in the form of Calcium Chloride, ( $\text{CaCl}_2$ ). This initiates several steps in the intrinsic pathway, leading to a fibrin clot. The time from addition of  $\text{CaCl}_2$  until clot formation is the APTT result in seconds.

## What Can Cause an Abnormal APTT?

One of the most important uses of an APTT reagent is to monitor the therapeutic coagulation inhibitors, like unfractionated heparin. Heparin assists antithrombin in its binding and inhibition of several coagulation factors. Antithrombin binds various active coagulation serine proteases (i.e., Factor Xa, Factor IIa [thrombin], Factor IXa) resulting in the slowing of clotting and prolongation of the APTT.

Different APTT reagents vary significantly in their sensitivity to unfractionated heparin.<sup>1</sup> A consensus has developed that the appropriate level of prophylactic heparin is 0.2 - 0.4 U/mL by protamine sulfate titration.<sup>2</sup> The appropriate response to these levels of heparin for APTT is considered to be between 1.5 - 2.5 times longer than the mid-point of the normal range.<sup>3,4</sup>

Significant deficiencies of procoagulant factors of the intrinsic pathway (XII, XI, IX, or VIII), those of the common pathway or inhibition of the reactions of these factors, can cause prolongation of the APTT. Reference ranges of the factors vary from approximately 50 - 150% of normal activity.<sup>5</sup>

Inhibitors affecting APTT include immunoglobulins that bind intrinsic or common pathway factors. One particular heterogenous group of immunoglobulins, collectively referred to as Lupus Anticoagulants, results in a variety of physiological manifestations and diagnostic indications. The APTTs of Lupus patients can vary considerably. Although the APTT is typically prolonged with Lupus, no APTT reagent has been shown to produce an abnormal result for all cases of Lupus.<sup>6</sup>

While deficiencies of fibrinogen and Vitamin K, as well as liver disease and oral anticoagulants, can cause prolonged APTT results, the situations in which the APTT is most importantly used are:

- Unfractionated Heparin Monitoring
- Screening for Coagulation factor deficiencies and their inhibitors
- Screening for Lupus Anticoagulants

## Raising the Standard for APTT Testing

IL provides a choice of solutions to satisfy an array of APTT testing requirements, from state-of-the-art synthetic micro-silica and ellagic acid technology to the classic silica formulation.





# Solutions

Improved Patient Care	Value	Ease of Use
<i>Provides consistent therapeutic ranges for improved heparin dosing and monitoring</i>	<i>Efficiently manages anticoagulated patients</i>	<i>Ensures precision of results and simplifies clinical interpretation through lot-to-lot consistency</i>
<i>Provides consistent sensitivity to intrinsic factor levels and Lupus Anticoagulants</i>	<i>Minimizes waste through reagent stability</i>	<i>Minimizes time-consuming parallel lot cross-over studies</i>
<i>Decreases risk of hemorrhagic and thrombotic complications by maintaining patients in therapeutic ranges</i>	<i>Standardizes performance for safer and more effective anticoagulant therapy</i>	<i>Liquid ready-to-use (except APTT S Lyo) and packaged with corresponding CaCl<sub>2</sub></i>



Instrumentation Laboratory is dedicated to providing unique solutions for today's modern hemostasis laboratories. We are committed to the highest quality products, new technologies through extensive ongoing research efforts, advanced system design and extensive service and support. The HemoSIL product line provides solutions for improved laboratory productivity and enhanced patient care.

\* Not available in all countries.

Providing the Right Solutions for Your APTT Testing

HemosIL Family of APTT reagents

APTT Reagent	APTT-SP	SynthASil	SynthAFax	APTT Lyo Silica***
Part Number and Size	0020006300 APTT reag. 5 x 9 mL 0.025M CaCl <sub>2</sub> 5 x 8 mL	0020006800 APTT reag. 5 x 10 mL 0.020M CaCl <sub>2</sub> 5 x 10 mL	0020007400 APTT reag. 5 x 10 mL 0.020M CaCl <sub>2</sub> 5 x 10 mL	0008468710 Cephalin 5 x 9 mL 0.025M CaCl <sub>2</sub> 5 x 8 mL
Liquid or Lyophilized	Liquid	Liquid	Liquid	Lyophilized
Activator	Silica	Colloidal Silica	Ellagic Acid	Micronized Silica
Phospholipids	Synthetic	Synthetic	Synthetic	Bovine brain cephalin
Reconstituted Stability at 2-8°C	APTT reagent 30 days in original vial CaCl <sub>2</sub> 30 days in original vial	APTT reagent 30 days in original vial CaCl <sub>2</sub> 30 days in original vial	APTT reagent 30 days in original vial CaCl <sub>2</sub> 30 days in original vial	APTT reagent 7 days in original vial CaCl <sub>2</sub> 30 days in original vial
Normal Range * **	24.4 - 35.0 seconds	25.4 - 38.4 seconds	19.7 - 27.6 seconds	24.9 - 36.8 seconds
Heparin Sensitivity	+++	+++	++	+++
Factor Sensitivity	+++	+++	++	++
Lupus Sensitivity	+++	++	+	+++

\* Due to many factors which can affect clotting times, each laboratory should establish its own reference range. \*\* Obtained on the ACL Family of instruments. \*\*\* Not available in the all countries.

Instrument Applications

	APTT-SP	SynthASil	SynthAFax	APTT Lyo Silica***
ACL TOP®	✓	✓	-	-
ACL Advance™	✓	✓	✓	✓
ACL ELITE/ELITE PRO - ACL 8/9/10000	✓	✓	✓	✓
ACL™ 100 - 7000	✓	✓	✓	✓

References

- Kitchen S, Preston F E. "The Therapeutic Range for Heparin Therapy: Relationship Between Six Activated Thromboplastin Time Reagents and Two Assays." *Thromb. Haemost.* 1996; 5: 734-9.
- Hirsh, J. "Heparin." *N. Eng. J. Med.* 1991; 324: 1565-74.
- Basu D, Gallus A, Hirsh J, Cade J. "A Prospective Study of the Value of Monitoring Heparin Treatment with the Activated \ Partial Thromboplastin Time." *N Engl J Med*, 1972; 287: 324-7.
- Hull, RD, Raskob GE, Pineo GF, Green D, Trowbridge AA, Elliott CG, Lerner RG, Hall J, Sparling T, Bretnell HR, Norton J, Carter CJ, Brant R, "Subcutaneous Low Molecular Weight Heparin Compared with Continuous Intravenous Heparin in the Treatment of Proximal-Vein Thrombosis." *N Engl J Med*, 1992; 326:975-82.
- Data on file, Instrumentation Laboratory.
- Triplet AD. "Laboratory Diagnosis of Lupus Anticoagulants." *Sem. Thromb.Haemost.* 1990; 16:183-192.

Werfen Group IVD

Worldwide Locations

Corporate Headquarters

Barcelona, Spain  
Tel. +34-93-4010101  
www.chwerfen.com

Instrumentation Laboratory Headquarters

Lexington, MA  
Tel. +1-781-861-0710  
www.ilus.com

US, Canada, Latin America

IL USA  
Lexington, MA  
Tel. +1-781-861-0710  
www.ilus.com

IL Canada  
Ontario, L4B 2N1  
Tel. +1-800-522-2025x6115

IL Mexico  
Col. Granada  
Tel. +52-55-5262-1760  
www.il-mexico.com.mx

Izasa Uruguay  
Montevideo  
Tel. +59-82-4818133

Pacific

Werfen Medical China  
Shanghai  
Tel. +86-21-32100745

Werfen Hong Kong  
Hong Kong  
Tel. +852-27927773

IL Japan  
Minato-ku, Tokyo  
Tel. +81-3-3437-6350

Werfen Medical Korea  
Seo Cho-ku  
Tel. +82-2-571-9246  
www.werfenmedical.com

Europe, Middle East, Africa

Werfen Austria  
Vienna  
Tel. +43-1-2565800-0

IL Belgium  
Zaventem  
Tel. +32-2-7252052  
www.il-be.com

Comesa Czech  
Prague  
Tel. +420-2-7816047

IL France  
Paris  
Tel. +33-1-53338600  
www.il-france.fr

IL Germany  
Munich  
Tel. +49-89-909070  
www.il-ger.de

IL Holland  
Breda  
Tel. +31-76-5480100  
www.il-nl.com

Comesa Hungary  
Budapest  
Tel. +36-1-4392910 or 11

IL Italy  
Milan  
Tel. +39-02-25221  
www.il-italia.it

Comesa Lithuania  
Kaunas  
Tel. +370-37-313157

Comesa Poland  
Warsaw  
Tel. +48-22-3361800

Izasa Portugal  
Camaxide  
Tel. +351-21-4247300  
www.izasa.com

Comesa Russia  
Moscow  
Tel. +7-495-9823723

Izasa Spain  
Barcelona  
Tel. +34-93-4010101  
www.izasa.com

IL UK  
Warrington, Cheshire  
Tel. +44-1925-81-0141  
www.il-uk.com

Instrumentation Laboratory, Comesa and Izasa are Companies of Werfen Group IVD

This Brochure is used throughout the world. Applications / tests listed may not yet be approved by the regulatory authorities in your country. Please contact the appropriate office above for specific information regarding availability in your country. IL product specifications are subject to modification to assure the highest quality performance. Some of the IL sites may still be in the process of completing ISO. ACL TOP, SynthASil, SynthAFax and HemosIL are registered trademarks and ACL is a trademark of Instrumentation Laboratory. © Instrumentation Laboratory 2003

