


# HemosIL<sup>®</sup> D-Dimer HS

## *DVT and PE Exclusion*

**Hemostasis**



**D-Dimer HS**

**D-Dimer Controls**

**Rapid and  
Cost-Effective for**  
*Enhanced Patient Care*

# D-Dimer in Clinical Practice

## In the Diagnosis of DVT and PE

Deep Vein Thrombosis (DVT) has an annual incidence of 0.5 - 1.2 per thousand. It is the primary cause of Pulmonary Embolism (PE), a potentially fatal event, occurring with an incidence of 0.2 - 0.6 per thousand. One diagnostic strategy for Venous Thromboembolism (VTE), which includes DVT and PE, begins with a clinical evaluation, followed by D-Dimer testing. Confirmatory tests generally involve imaging techniques. For DVT, these are lower limb venous compression ultrasonography (CUS) and venography, which is invasive and considered the gold standard. For PE, imaging tests are helical computerized tomography (CT) scanning or ventilation/perfusion (V/Q) lung scanning and angiography, which is also invasive.

A substantial number of publications in the past few years report the use of D-Dimer, together with pre-test probability (PTP) assessment, as a safe, cost-effective management strategy for the evaluation of patients presenting to emergency departments with clinically suspected VTE. This approach allows DVT and/or PE to be ruled out in outpatients with suspected VTE with low or low-moderate PTP and a negative D-Dimer, hence reducing the number of imaging tests required, particularly those invasive.

## In the Diagnosis of DIC


Disseminated Intravascular Coagulation (DIC) is considered a systemic thrombohemorrhagic disorder. It is a secondary manifestation to an underlying disorder, such as sepsis, trauma or malignancy. DIC has been described as low-grade, compensated or fulminant, but when a patient presents with suspected DIC, such a clear differentiation of the disease is not apparent. The patient may be at any stage between these extremes with progressing severity. Therefore, early diagnosis and appropriate treatment are of primary importance in DIC diagnosis. Recent definitions of clinical laboratory criteria for the diagnosis and staging of DIC include D-Dimer measurement as an important contributor to the algorithm for the DIC scoring system.

## Other Clinical Applications

D-Dimer testing has been evaluated in a number of clinical applications such as predictive factor for recurrences of VTE after discontinuation of oral anticoagulant therapy, and as an indicator of pregnancy complications. Other studies suggest that D-Dimer levels in the normal population may indicate the risk for arterial thrombosis.



**D-Dimer Levels in Patient Groups and Normal Subjects**



An example of D-Dimer levels in different clinical conditions: normal patients, PE, DVT, DIC.




## HemosIL D-Dimer HS

HemosIL D-Dimer HS is a fully automated, D-Dimer assay for the exclusion of DVT and PE, in conjunction with a clinical Pretest Probability (PTP) assessment model.

The latex reagent is a suspension of polystyrene latex particles of uniform size, coated with the F(ab')<sub>2</sub> fragment of a monoclonal antibody, highly specific for the D-Dimer domain included in fibrin soluble derivatives. The use of the F(ab')<sub>2</sub> fragment allows for a more specific D-Dimer detection, avoiding the interference from some endogenous factors like the Rheumatoid Factor.

The assay is optimized for a wavelength of 671 nm, providing very high analytical sensitivity, enhanced working range and minimized interferences from hemoglobin.

In multi-center management studies with over 800 samples from patients suspected of VTE, HemosIL D-Dimer HS demonstrated 100% Negative Predictive Value (NPV) and 100% Sensitivity at a cut-off value of 230 ng/mL.



## Analytical Characteristics on the ACL TOP® Hemostasis Testing System

The data were obtained according to Clinical Laboratory Standards Institute (CLSI) guidelines.

<b>Linearity</b>		
<i>Up to 3,680 ng/mL without rerun</i>		
<i>Up to 69,000 ng/mL with automated rerun</i>		
<b>Imprecision</b>		
	<b>Mean (ng/mL)</b>	<b>CV (Total)</b>
<i>Low D-Dimer Control</i>	314	7%
<i>High D-Dimer Control</i>	677	7%
<b>Detection Limit</b>		
<i>21 ng/mL</i>		
<b>Interferences</b>		
<i>Hemoglobin</i>	<i>none up to 500 mg/dL</i>	
<i>Bilirubin</i>	<i>none up to 18 mg/dL</i>	
<i>Triglycerides</i>	<i>none up to 1,327 mg/dL</i>	
<i>Rheumatoid Factor</i>	<i>none up to 1,400 IU/mL</i>	
<b>Stability of Opened / Reconstituted Reagents</b>		
<i>Latex Reagent</i>	<i>1 month at 2-8°C</i>	
<i>Reaction Buffer</i>	<i>1 month at 2-8°C</i>	
<i>Calibrator</i>	<i>1 month at 2-8°C</i>	



# VTE Management Studies with HemosIL D-Dimer HS

## Enhancing Diagnosis of DVT and PE

Three hundred frozen samples from patients admitted consecutively to the Emergency Department with suspected PE or DVT (frequency of venous thromboembolism 26%) were evaluated during the single center study. Seventy-eight were confirmed as VTE-positive (47 PE and 31 DVT) by standard objective tests and the remaining 222 were

confirmed as negative. The optimal cut-off value on the ACL TOP Hemostasis Testing System was determined to be 230 ng/mL. Patient samples were also tested with a traditional automated, latex-based method and an ELISA assay. The cut-off values were those recommended by the manufacturer. The results are summarized in the table below.

	<i>HemosIL D-Dimer HS</i>	<i>Traditional Automated Latex</i>	<i>ELISA Method</i>
<i>Sensitivity (95% CI)</i>	100% (95.4 - 100%)	100% (95.1 - 100%)	100% (95.4 - 100%)
<i>Specificity (95% CI)</i>	46.8% (40.1 - 53.6%)	35.9% (29.6 - 42.6%)	34.7% (28.4 - 41.3%)
<i>NPV (95% CI)</i>	100% (96.5 - 100%)	100% (95.4 - 100%)	100% (95.3 - 100%)

## Fewer False Positive Results Versus Traditional Latex and ELISA Methods

The management study comparing HemosIL D-Dimer HS to a traditional automated, latex-based assay and an ELISA

assay demonstrated the following results in the sample population tested.

### Distribution of Patient Samples from the Management Study

	<i>HemosIL D-Dimer HS</i>	<i>Traditional Automated Latex</i>	<i>ELISA Method</i>
<i>Total number of patients with suspected VTE, tested with D-Dimer</i>	300	294*	300
<i>Number of patients excluded for VTE (D-Dimer &lt; cut-off)</i>	104	79	77
<i>Number of patients undergoing investigation with imaging techniques (D-Dimer &gt; cut-off)</i>	196	215	223
<i>Outcome:</i>			
<i>Confirmed for DVT</i>	31	31	31
<i>Confirmed for PE</i>	47	43*	47
<i>False Negative for DVT or PE</i>	0	0	0
<i>False Positives for VTE</i>	118	141	145

\* Six samples did not report any result due to instrument errors and were excluded from the analysis; four were PE and two were normal samples.

HemosIL D-Dimer HS demonstrated fewer false positive patient samples (N=118) than the traditional automated, latex-based assay (N=141) and the ELISA method (N=145);

thus, supporting the value of the assay as a cost-effective tool for the management of VTE patients.

# Excluding DVT and PE



## Multi-center Studies

To confirm the 230 ng/mL cut-off value, a multi-center study was performed at four hospitals on 668 samples from patients admitted consecutively to the Emergency Department with suspected PE or DVT. As part of the study, patients underwent a PTP assessment using the Wells model and were classified as having a high, moderate, or low probability of PE or DVT. Patients with a negative D-Dimer test result and a low PTP score underwent no further diagnostic testing and were followed-up after three months for development of PE or DVT. For patients with a negative

D-Dimer test result and a moderate PTP, the physicians determined if a three-month follow-up or imaging techniques were required. Patients with a positive D-Dimer test result or a high PTP score underwent imaging techniques. The overall prevalence in the total population of samples was 16.1% (58/361) for PE and 20.2% (62/307) for DVT. The sensitivity, specificity and NPV of HemosIL D-Dimer HS, using the previously established clinical cut-off of 230 ng/mL, is summarized below for DVT and PE with the corresponding 95% confidence intervals (CI).

DVT Population	All samples	High PTP	Low + Moderate PTP
<i>n</i>	307	54	253
Sensitivity (95% CI)	100% (62/62) (94.2 - 100%)	100% (28/28) (87.7 - 100%)	100% (62/62) (89.7 - 100%)
Specificity (95% CI)	38.4% (94/245) (32.2 - 44.8%)	34.6% (9/26) (17.2 - 55.7%)	38.8% (62/62) (32.3 - 45.6%)
NPV (95% CI)	100% (94/94) (96.2 - 100%)	100% (9/9) (66.4 - 100%)	100% (62/62) (95.8 - 100%)

PE Population	All samples	High PTP	Low + Moderate PTP
<i>n</i>	361	28	333
Sensitivity (95% CI)	100% (58/58) (93.8 - 100%)	100% (10/10) (62.2 - 100%)	100% (48/48) (92.6 - 100%)
Specificity (95% CI)	35.6% (108/303) (30.2 - 41.3%)	16.7% (3/18) (3.6 - 41.4%)	36.8% (105/285) (31.2 - 42.7%)
NPV (95% CI)	100% (108/108) (96.6 - 100%)	100% (3/3) (29.9 - 100%)	100% (105/105) (96.5 - 100%)

## HemosIL D-Dimer HS Safely Excludes DVT and PE when Used in Conjunction with a PTP Assessment Model

Using a 230 ng/mL cut-off, HemosIL D-Dimer HS successfully excluded DVT and PE in both studies' patient populations with an NPV of 100% and sensitivity of 100%.

### Single Center Management Study

- 300 outpatients
- NPV = 100%
- Specificity = 46.8%

### Multi-Center Management Study

- 668 outpatients
- NPV = 100%
- Specificity for DVT = 38.4%
- Specificity for PE = 35.6%

# HemosIL D-Dimer HS

## DVT and PE Exclusion

<b>D-Dimer HS Kit Composition P/N 0020007700</b>	<b>3 x 2 mL Latex Reagent (Lyo)</b>
	<b>3 x 8 mL Reaction Buffer (Liq)</b>
	<b>2 x 1 mL Calibrator (Lyo)</b>
<b>D-Dimer Controls P/N 0020008610</b>	<b>5 x 1 mL Low D-Dimer Control</b>
	<b>5 x 1 mL High D-Dimer Control</b>



### References

- Perrier A, Bounameaux H. Cost-effective diagnosis of deep vein thrombosis and pulmonary embolism. *Thromb Haemost* 2001; 86:475-487.
- Anderson DR, Wells PS. Improvement in the diagnostic approach for patients with suspected deep vein thrombosis or pulmonary embolism. *Thromb Haemost* 1999; 82(2): 878-886.
- Wells PS, Anderson DR, Rodger M, Forgie M, Kearon C, Dreyer J, Kovacs G, Mitchell M, Lewandowski B, Kvacs MJ. Evaluation of D-Dimer in the diagnosis of suspected deep-vein thrombosis. *N Engl J Med* 2003; 349: 1227-1235.
- Bockensted P. D-Dimer in venous thromboembolism. *N Engl J Med* 2003; 349: 1203-1204.
- Oswal CT, Menon V, Stouffer GA. The use of D-Dimer in emergency room patients with suspected deep vein thrombosis: a test whose time has come. *J Thromb Haemost* 2003; 1: 635-636.
- Anderson DR, Kovacs MJ, Kovacs G, Stiell I, Mitchell M, Khoury V, Dryier J, Ward J, Wells PS. Combined use of clinical assessment and D-Dimer to improve the managements of patients presenting to the emergency department with suspected deep vein thrombosis (the EDITED study). *J Thromb Haemost* 2003; 1: 645-651.
- Curtin N, Highe G, Harris M, Braunstein A, Demattia F, Coss L. Extensive evaluation of the Instrumentation Laboratory IL Test D-Dimer immunoturbidimetric assay on the ACL 9000 determined the D-Dimer cut-off value for reliable exclusion of venous thromboembolism. *Laboratory Hematology* 2004; 10:88-94.
- Legnani C, Pancani C, Palareti G, Guazzaloca G, Coccheri S. Performance of a new, fast D-Dimer test (IL Test D-Dimer) for the management of outpatients with suspected deep vein thrombosis in emergency situations. *Fibr & Proteol* 1999; 13(3): 139-141.
- Kovacs MJ, MacKinnon KM, Andreson D, O'Rourke K, Keeny M, Kea C, Ginsberg J, Wells PS. A comparison of three rapid D-Dimer methods for the diagnosis of venous thromboembolism. *Br J Haematol* 2001; 115(1):140-144.
- Taylor FB Jr, Toh CH, Hoots WK, Wada H, Levi M. Toward definition, clinical and laboratory criteria, and scoring system for disseminated intravascular coagulation. *Thromb Haemost* 2001; 86: 1327-1330.
- Hambleton J, Leung LL, Levi M. Coagulation: Consultative hemostasis. *Hematology (Am Soc Hematol Educ Program)* 2002: 335-352.
- Morse M. Establishing a normal range for D-Dimer levels through pregnancy to aid in the diagnosis of pulmonary embolism and deep vein thrombosis. *J Thromb Haemost* 2004; 2: 1202-1204.
- Palareti G, Legnani C, Cosmi B, Valdre L, Lunghi B, Bernardi F, Coccheri S. Predictive value of D-dimer test for recurrent venous thromboembolism after anticoagulation withdrawal in subjects with a previous idiopathic event and in carriers of congenital thrombophilia. *Circulation* 2003; 108(3): 313-318.
- Nieuwenhuizen W. A Reference Material for Harmonisation of D-Dimer Assays. *Thromb Haemost*. 1997; 77 (5): 1031-1033.
- de Moerloose P, Vanrusselt M, Reber G and Arnout J. Performances of the HemosIL D-Dimer HS assay for the exclusion of venous thromboembolism. *JThromb Haemost* 2005; 3, 2361-3
- Legnani C, Cosmi B, Cini M, Mattarozzi S, Lo Manto G, Poggi M, Palareti G. Effects of different cut-off values of three quantitative D-Dimer methods to establish duration of oral anticoagulation treatment after a first idiopathic episode of venous thromboembolism, XIX SISET Congress, Milan, Italy, September 2006. Poster session; P152.

## Werfen Group IVD

### Worldwide Locations

#### Corporate Headquarters

Barcelona, Spain  
Tel. +34-93-4010101  
www.werfengroup.com

#### Instrumentation Laboratory Headquarters

Lexington, MA  
Tel. +1-781-861-0710  
www.ilwww.com

#### US, Canada, Latin and South America

**IL USA**  
Lexington, MA  
Tel. +1-781-861-0710  
www.ilus.com

**Werfen Medical Brazil**  
São Paulo, SP  
Tel. +55-11-37353641

**IL Canada**  
Richmond Hill, ON  
Tel. +1-800-522-2025 x6115

**IL Mexico**  
Col. Granada  
Tel. +52-55-5262-1760  
www.il-mexico.com.mx

**Izasa Uruguay**  
Montevideo  
Tel. +59-82-4818133

#### Pacific

**Werfen Medical China**  
Shanghai  
Tel. +86-21-32100745

**Werfen Hong Kong**  
Hong Kong  
Tel. +852-27927773

**IL India**  
Janakpuri, New Delhi  
Tel. +91-11-25510137

**IL Japan**  
Minato-ku, Tokyo  
Tel. +81-3-3437-6350

**Werfen Medical Korea**  
Seo Cho-ku  
Tel. +82-2-571-9246  
www.werfenmedical.com

#### Europe, Middle East, Africa

**Werfen Austria**  
Vienna  
Tel. +43-1-2565800-0

**IL Belgium**  
Zaventem  
Tel. +32-2-7252052  
www.il-be.com

**Comesa Czech**  
Prague  
Tel. +420-2-7816047

**IL France**  
Paris  
Tel. +33-1-53338600  
www.il-france.fr

**IL Germany**  
Munich  
Tel. +49-89-909070  
www.il-ger.de

**IL Holland**  
Breda  
Tel. +31-76-5480100  
www.il-nl.com

**Comesa Hungary**  
Budapest  
Tel. +36-1-4392910 or 11

**IL Italy**  
Milan  
Tel. +39-02-25221  
www.il-italia.it

**IL Lithuania**  
Kaunas  
Tel. +370-37-313157

**Comesa Poland**  
Warsaw  
Tel. +48-22-3361800

**Izasa Portugal**  
Camaxide  
Tel. +351-21-4247300  
www.izasa.com

**IL Russia**  
Moscow  
Tel. +7-495-9823723

**Izasa Spain**  
Barcelona  
Tel. +34-93-4010101  
www.izasa.com

**IL UK**  
Warrington, Cheshire  
Tel. +44-1925-81-0141  
www.il-uk.com

Instrumentation Laboratory, Comesa and Izasa are Companies of Werfen Group IVD

Applications/tests listed may not yet be approved by the regulatory authorities in your country. Please contact the appropriate office above for specific information regarding availability in your country. IL product specifications are subject to modification to assure the highest quality performance. Some of the IL sites may still be in the process of completing ISO. ACL TOP and HemosIL are registered trademarks of Instrumentation Laboratory.  
© Instrumentation Laboratory 2005

